

EMERGENCY BAGS AND RUCKSACKS

Boscarol's bags and rucksacks are the result of years of study and cooperation with doctors and rescuers, with experts in the field of sea and alpine rescue teams besides international civil protection. Opposite to what is often thought, our medical bags and rucksacks are specially manufactured through particular processes that include the choice of the material which must be resistant to mechanical strain, specific sewing and block-free zippers easy to open although self protection devices are involved. Nothing is to be considered casual, not even the choice of the colour and the fantasy of the fabric used for the pediatric line.

Boscarol company is also proud to be the manufacturer of the original RESCUE BAG, designed and produced since 1985 by the founder of the company itself Mr. Oscar Boscarol, which is currently the most copied bag by the competitors worldwide.

The wide range of bags available aims to satisfy all of our customers needs, thus helping to make our brand worldwide famous.



Emergency rucksacks

Tear-proof impermeable fabric material
Double handle for upper or lateral carrying
Wide frontal pocket total opening with 2 internal elastics strips for equipment fixing
3 perimeter external pockets
Lateral and frontal refractive strips
Supplied with 8 coloured removable pouches with transparent frontal and Velcro-strips
Central bulkhead with zip, complete with elastics
Reinforced bottom basin and rubber feet
Dimensions: about 47x27x48h cm
Weight: about 3 kg



BOR87000 RESCUE RUCKSACK VIPER

Emergency rucksack MINI Viper
Made of red tear-proof impermeable fabric material
Double handle for upper or lateral carrying
Wide frontal pocket with 2 lines of elastics suitable for equipment fixing inside the pocket
3 perimeter external pockets
Lateral, frontal and backside refractive strips
Featured with 4 coloured removable pouches with transparent frontal and Velcro-strips
Reinforced bottom basin and rubber feet
Dimensions: about 43x23x47h cm
Weight: about 2,7 kg
Supplied empty



BOR87002 RESCUE RUCKSACK MINI VIPER

Emergency rucksack Viper in Tyrpol
Made of red tear-proof impermeable fabric material
Equipped with upper handle and shoulder straps for transport
Wide frontal total-opening pocket with 2 internal elastic strips for equipment fixing
2 perimeter external pockets
Lateral and frontal yellow refractive strips
Supplied with 6 coloured removable pouches with transparent frontal window and Velcro-strips
Central bulkhead with zip, complete with elastics
Reinforced bottom basin and rubber feet
Dimensions: 34,5x19,5x46,5 cm
Supplied empty



BOR87012 RESCUE RUCKSACK VIPER TYRPOL



Tearproof impermeable red Tyrol material
2 frontal equipped pockets
2 empty external perimeter pocket 1 of which with
4 coloured pouches with transparent frontal
Inner rigid, shockproof basin with oxygen holder
4 inner removable coloured pouches with transparent frontal
Dimensions: about 48x35x54 h cm – 54lts
Weight: about 4,6 kg

BOR86030

RESCUE RUCKSACK MAXI 02 TYRPOL RED
EMPTY



Tearproof impermeable material
2 frontal pockets and 2 yellow refractive strips
8 internal coloured pouches with transparent frontal
4 of which with transport handle
Central oxygen cylinder holder straps
Internal perimeter provided with elastics
Dimensions: about 35x30x53h cm
Weight 4,5 kg

BOR87020

RUCKSACK MASTERMED EMPTY



Tearproof impermeable material
Frontal pocket for association logo
2 external perimetric pockets
Book opening
Inner elastics and pockets for equipment fixing
Rigid bottom side
Dimensions: about 45x24x48(h) cm
Weight: about 2,4 kg

BOR80200/P

PC RUCKSACK EMPTY



Emergency rucksack Maximed empty
Made of red tearproof impermeable fabric material
Padded straps
2 frontal equipped pockets
1 perimeter external pocket
1 perimeter external pocket with 4 coloured removable pouches
4 inner removable coloured pouches with transparent frontal
Elastics alongside the internal perimeter of the rucksack
Oxygen cylinder holder straps placed centrally
Dimensions: about 48x30x52(h) cm – 54 lts
Weight: about 4,5 kg

BOR87010

ZAINO MAXIMED VUOTO

ATLS rucksacks



Tearproof impermeable material
Padded shoulder straps
1 frontal transparent pocket for association logo
1 wide frontal pocket
Inner, rigid basin with compartments
8 coloured removable pouches with transparent frontal, 4 of which with handle
Dimensions: about 32x25x53(h) cm – 35 lts
Weight: about 4 kg

BOR81146 ATLS MODULE RUCKSACK EMPTY



Tearproof impermeable material
Padded shoulder straps
1 frontal transparent pocket for association logo
1 wide frontal pocket
Inner, rigid basin with central oxygen bottle holder
4 coloured removable pouches with transparent frontal
Dimensions: about 32x25x53(h) cm – 35 lts
Weight: about 3,3 kg

BOR81145 EMERGENCY RUCKSACK ATLS O2 EMPTY

BOR81147 EMERGENCY RUCKSACK ATLS WK EMPTY

Emergency Trolley

Trolley Rucksack Master empty
2 reflective strips
2 frontal equipped pockets
Inner rigid, shockproof basin with central support for oxygen cylinder
4 removable coloured pouches with transparent frontal
Dimensions: about 33x35x59(h) cm – 43 lts
Weight: 5 kg

BOR87220 RUCKSACK MASTER TROLLEY TYRPOL EMPTY



TRAUMA bags



Tearproof and impermeable red fabric material
4 reflective straps and totally removable cover
2 external lateral pockets
2 external lateral pockets with elastics and additional pockets
4 inner modular spaces
1 inner removable pouch with transparent frontal
Dimensions: about 45x33x30(h) cm
Weight: about 2,2 kg

BOR89040 TRAUMA BAG 10 RED EMPTY



Tearproof and impermeable red fabric material
4 reflective straps and totally removable cover
2 external lateral pockets
2 external lateral pockets with elastics and additional pockets
5 inner modular spaces
1 inner removable pouch with transparent frontal
Dimensions: about 53x33x34(h) cm
Weight: about 2,2 kg

BOR89050 TRAUMA BAG 20 RED EMPTY

Functional, capacious, resistant and compact emergency bag
Made of tearproof and impermeable red fabric material
4 reflective straps
Totally removable cover
Easy to carry either by hand, on shoulder or as a backpack
2 external lateral pockets
3 external lateral pockets with elastics and additional pockets
5 inner modular spaces
2 inner removable pouches with transparent frontal
Dimensions: about 63x37x31(h) cm
Weight: about 2,3 kg



BOR89060 TRAUMA BAG 30 RED EMPTY

Emergency bags



Tearproof and impermeable fabric material
1 capacious outer pocket
5 inner adjustable spaces
Inside provided with elastics and a documents pocket
1 transparent frontal pocket for association logo
Dimensions: about 36x18x35(h) cm
Weight: about 1,6 kg

BOR87300 BAG 33 EMPTY



Impermeable tearproof material
To be carried by handle and shoulder strap
3 transparent pockets inside inner top
Inner rigid basin
Dimensions: about 48x26x16h cm
Weight: about 1,2 kg

BOR80390 MULTICASE RED EMPTY



Tearproof and impermeable material
5 dividers with pocket
Belt loop of about 5,5 cm
Dimensions: about 14x3x19(h) cm
Complete with:
1 first aid scissors
1 anatomic forceps
Cocher forceps
1 Pen Light
1 Tourniquet

BOR84064 E-HOLSTER BLUE COMPLETE

BOR84065 E-HOLSTER RED COMPLETE

BOR84062 E-HOLSTER BLUE EMPTY

BOR84063 E-HOLSTER RED EMPTY



Particularly functional belt case
Made in blue tearproof impermeable fabric material
Suitable to contain one sphygmomanometer
3 inner spaces with closing
Dimensions: about 34x11x8(h) cm
Weight: about 100 gm

BOR205 HOLSTER PROFESSIONAL ROSSO

BOR206 HOLSTER PROFESSIONAL BLU

RESCUE and PROFI bags



Tyrol material
2 frontal transparent pockets for association logo
4 inner coloured pouches
Upper part with inner elastics suitable for equipment fixing
Bottom side with shockproof inner basin
Arranged to feature one oxygen cylinder
Dimensions: about 48x26x28(h) cm
Weight: about 2,2 kg

BOR80340 RESCUE BAG MODULE



Tearproof impermeable material
2 frontal transparent pockets for any association logo
Bottom side with inner shockproof basin
Upper side with inner fixed elastics for material fixing
Dimensions: about 39x18x25(h) cm
Weight: about 1 kg

BOR80290 MINI RESCUE BAG RED EMPTY



Impermeable tearproof material
2 frontal transparent pockets for association logo
Bottom side with shockproof inner basin
Arranged to feature one oxygen cylinder
Upper part with inner elastics suitable for equipment fixing
Dimensions: about 48x26x28(h) cm
Weight: about 1,8 g

BOR80310 RESCUE BAG RED INTERNAL OPENING EMPTY



Impermeable and tearproof red fabric material
Modular interior
1 frontal pocket
Arranged to hold 1-2 lts oxygen cylinder
2 wide lateral external pockets and internal pocket
Dimensions: about 60x25x27(h) cm
Weight: 1,4 kg
PROFI BAG MINI
2 wide lateral external pockets
Internal pocket with transparent frontal
Dimensions: about 36x26x20(h) cm
Weight: 0,7 kg

BOR89010 PROFI BAG EMPTY

BOR89000 PROFI BAG MINI EMPTY

Pediatric bags



Tearproof and water-resistant material
Shockproof rigid internal basin with 2 compartments
2 removable pouches with transparent frontal
On top, 3 inner transparent pockets
Dimensions: about 40x18x13(h) cm
Weight: about 0,9 kg

BOR89100 BABY BAG EMPTY



Tearproof and water-resistant material
Completely removable cover
3 lateral external pockets
1 lateral external pocket with elastics
4 inner modular spaces
1 removable pouch with transparent frontal
Dimensions: about 45x33x30(h) cm
Weight: about 1,8 kg

BOR89090 BABY CHILDREN REANIBAG EMPTY



Tearproof impermeable Tyropol material
Shoulder straps with sectorially-coloured measurement system for easy identification
Divided into 3 compartments
8 coloured removable pouches with transparent frontal
Dimensions: about 40x24x37h cm
Weight: about 2,7 Kg

BOR86064 BABY TRIS COLOURED EMPTY

Oxygen bags



Padded oxygen cylinder bag (up to 5 litres)
Tearproof impermeable material
External pockets for accessories
Equipped with handle and shoulder strap
Dimensions: about 75x35x17 cm
Weight: about 0,5 kg

BOR80614 OXYGEN CYLINDER BAG 5LT EMPTY



Padded oxygen cylinder holder
High-resistance impermeable TYRPOL material
Equipped with handle and stretcher safe fixing system
Arranged to hold both 1 or 2 lts oxygen cylinders
Dimensions: about 55x25x14 cm
Weight: about 0,9 kg

BOR80612 TYRPOL OXIBAG EMPTY



Oxygen cylinder bag (up to 2 litres)
Nylon impermeable tearproof material
Dimensions: Ø 10x36 cm

BOR80615 OXYGEN CYLINDER BAG 2 lt EMPTY



Padded oxygen cylinder bag (up to 2 lts)
Tearproof impermeable material
Equipped with handle and shoulder strap
Equipped with compartments, pockets and elastics for items fixing
Dimensions: about 51x22 cm
Weight: about 0,7 kg

BOR80613 OXYGEN CYLINDER BAG 2LT EMPTY

Pouches



DIM. 25x25x8 cm

BOR81042	DOUBLE BLUE POUCH
BOR81044	DOUBLE RED POUCH WITH ELASTICS



DIM. 31x15x8 cm

BOR89450	BLUE BIG STANDARD POUCH
BOR89452	RED BIG STANDARD POUCH
BOR89454	GREEN BIG STANDARD POUCH
BOR89456	YELLOW BIG STANDARD POUCH



DIM. 28x12x9(h) cm
 DIM. 31x15x7(h) cm
 DIM. 25x25x8(h) cm
 DIM. 25x25x8(h) cm

BOR81050	SMALL FANTASY POUCH
BOR81056	DOUBLE FANTASY POUCH
BOR81052	FANTASY POUCH MEDIUM
BOR81054	FANTASY POUCH DOUBLE WITH ELASTIC



DIM. 35x10x13 cm

BOR89468	BLACK BIG POUCH
----------	-----------------



DIM. 19x12x6 cm

BOR81058	VIPER YELLOW POUCH
BOR81060	VIPER GREEN POUCH
BOR81062	VIPER RED POUCH
BOR81064	VIPER BLUE POUCH
BOR81066	VIPER VIOLET POUCH
BOR81068	VIPER BLACK POUCH
BOR81070	VIPER ORANGE POUCH



DIM. 30x12x13 cm

BOR89440	GREEN BIG POUCH WITH HANDLE
BOR89442	RED BIG POUCH WITH HANDLE
BOR89444	BLUE BIG POUCH WITH HANDLE
BOR89446	YELLOW BIG POUCH WITH HANDLE



DIM. 27x12x9 cm

BOR81024	STANDARD PINK POUCH
BOR81028	STANDARD RED POUCH
BOR81030	STANDARD GREEN POUCH
BOR81032	STANDARD YELLOW POUCH
BOR81034	STANDARD BLUE POUCH
BOR81036	STANDARD VIOLET POUCH
BOR81038	STANDARD ORANGE POUCH
BOR81040	STANDARD BLACK POUCH



DIM. 23x9x8 cm

BOR81082	BABY TRIS POUCH RED
BOR81084	BABY TRIS POUCH GREEN
BOR81086	BABY TRIS POUCH YELLOW
BOR81088	BABY TRIS POUCH BLUE
BOR81089	BABY TRIS POUCH NAVY BLUE
BOR81090	BABY TRIS POUCH PURPLE
BOR81092	BABY TRIS POUCH ORANGE
BOR81094	BABY TRIS POUCH BLACK



Made in tearproof impermeable red material
Book opening
Complete with 3 net pockets for emergency equipment
Waist wearing mode
Dimensions: about 13x6x18(h) cm
Weight: about 300 gm

BOR220	FIRST AID "210" VUOTO
--------	-----------------------



MIL1370	BIG POUCH OLIVE GREEN
MIL1372	SMALL POUCH OLIVE GREEN



Empty coloured mini pouch with transparent frontal
Tearproof impermeable black fabric
Velcro strips for fixing inside bags and rucksacks
Dimensions: about 15x12x8,5(h) cm

BOR89458	BLACK MINI POUCH
BOR89460	RED MINI POUCH
BOR89462	BLUE MINI POUCH
BOR89464	GREEN MINI POUCH



Fabric empty bag
Tearproof impermeable blue material
With shoulder strap
Dimensions: about 29x16x16 cm

BOR84072	FABRIC BAG
----------	------------



Empty coloured pouch with transparent frontal
Tearproof impermeable fabric
With handle
Velcro strips for fixing inside bags and rucksacks
Suitable for Rescue Bag
Dimensions: about 23x14x11(h) cm

BOR89420	YELLOW POUCH WITH HANDLE
BOR89422	BLUE POUCH WITH HANDLE
BOR89424	GREEN POUCH WITH HANDLE
BOR89426	RED POUCH WITH HANDLE

Pouches and aluminium cases



Small bag made in PVC
Washable
Tearproof impermeable red fabric
Designed to be attached to the OB2012 and OB1000
Suction Unit transport bags
Suitable to contain small sets
Dimensions: 10x7x13(h) cm

BOR89469 MINI PVC RED BAG

BOR89471 MINI PVC BLUE BAG



Emergency medical case empty in aluminium
Made of anodized aluminium
High resistance and long durability
Inside made in washable and disinfectable thermo-format with modular system
Arranged to accommodate oxygen cylinders (1 or 2 litres)
Type B: removable internal basin made in thermo-format with rail ampoules holder (max 60 pcs). Dimensions: about 52x36x24(h) cm. Weight: about 7,6 kg
Type A: Dimensions: about 52x36x17(h) cm
Weight: about 5,7 kg

VAL42000 EMERGENCY MEDICAL CASE ALUMINIUM TYPE A

VAL42010 2 AMPOULE HOLDER RAIL DEPOT FOR CASE TYPE A

VAL43000 EMERGENCY MEDICAL CASE ALUMINIUM TYPE B



Carbon steel cabinet, coated in epoxy polyester
Rounded edges
Three steel shelves
Door with 3 accessories holder
Locker with key
Weight: 6,60
Dimensions: 402x202x483(h) mm

CAV605 EMPTY CABINET



Metal Cabinet for Defibrillator with see through door
Made of Galvanized Carbon Steel PO2
Enameled by Epoxy Polyester white powder RAL9010 Heated at 180°
Transparent door with opening at 180° with a galvanized strong pivot
Strong Nickel plated zinc alloy closing device. Can be sealed or locked
Back holes for wall hanging.
Weight empty item: 6 kg
Size: 425x425x160 mm

DEF040 VISIO DEFIBRILLATOR METAL CABINET

Emergency cases



Rigid Case Explorer
Water repellent up to 10m depth
Pressurization valve
Resistant to corrosive agents
Strong click closing predisposed for padlock or lead

VAL43151	EXPLORER CASE MINI - DIM: 305x195x270(h) mm
VAL43152	EXPLORER CASE P - DIM: 410x205x340(h) mm
VAL43162	EXPLORER CASE M - DIM: 520x230x435(h) mm
VAL43164	EXPLORER CASE EMPTY - DIM: 445x345x190(h) mm
VAL43172	EXPLORER CASE G - DIM: 650x240x510(h) mm
VAL43174	CASE EXPLORER P WITH 4 EMPTY POUCHES DIMENSIONI 410x205x340(h) mm
VAL43176	CASE EXPLORER M WITH 5 EMPTY POUCHES DIMENSIONI 520x230x435(h) mm
VAL43178	CASE EXPLORER G WITH 6 EMPTY POUCHES DIMENSIONI 650x240x510(h) mm



San type case empty
Made in orange shockproof ABS material
Wall bracket anchorage system
2 equal sides made of 2 adjustable
inner compartments
2 transparent dividers
Dimensions: about 32,5x12,5x21,5(h) cm
Weight: about 1,1 Kg

VAL40130 SAN CASE IN ABS EMPTY



Multi type case empty
Made in orange shockproof ABS material
Wall bracket anchorage system
2 equal sides made of 2 adjustable inner compartments
2 transparent dividers
Dimensions: about 42x15x30(h) cm
Weight: about 3,3 Kg

VAL40110 MULTI CASE IN ABS EMPTY

Boscarol High Impact Case

Boscarol High Impact orange case with high resistant thick walls
 Automatic pressure release valve
 Easy to open double throw latches
 Load tested soft grip handle
 Holes for double padlocking
 Holes for shoulder strap
 Stackable interlocking system
 Thick body material
 Extra fixing points
 Fully length hinges
 Nylon pins, corrosion free
 Pick N Pluck foams in the bottom and convoluted foam available on demand
 Screen printing space on the lid
 Rectangular resin label space on the lid
 Weight: 1,75 kg



VAL43140	BOSCAROL HIGH IMPACT CASE P - 33,6x30x14,8(h) cm
VAL43141	BOSCAROL HIGH IMPACT CASE M - 46,4x36,6x17,6(h) cm
VAL43142	BOSCAROL HIGH IMPACT CASE G - 55,5x42,8x21,1(h) cm



Shockproof orange ABS material
 Equipped with wall bracketing
 Dimensions: about 39,5x13,5x26(h) cm
 Weight: about 1,6 kg

CAV401	MEDIC 2 CASE ORANGE COVER
CAV399	MEDIC 2 CASE TRANSPARENT COVER



Flexi case empty
 Made in shockproof orange polypropylene material
 Equipped with wall bracketing
 Dimensions: about 25,5x9x19,5 (h) cm
 Weight: about 0,45 kg

CAV406	FLEXI CASE EMPTY
--------	------------------

Drugs bags

Book-opening drugs case
Tearproof impermeable red material
2 removable red plates suitable for 60 ampoules
Interior padded dividers with pockets, for needles and further medical material
2 Velcro-strips for the fixing in bags and emergency rucksacks
Dimensions: about 32x22,5x7,5h cm
Weight: 1,2 kg

BOR89310	RAIL DRUGS CASE
BOR89311	EASY RAIL DRUGS CASE



Book-opening drugs case empty
Made in tearproof impermeable red fabric material
2 Velcro-strips for the fixing in bags and emergency rucksacks
2 removable red plates with adjustable elastics suitable for 72 ampoules
Padded inner divider with pockets for needles and additional medical equipment
Dimensions: about 32x9x23(h) cm
Weight: about 0,8 kg

BOR89316	DRUGS CASE WITH ELASTICS
BOR89326	MAXI DRUGS CASE WITH ELASTICS 45X9,5X33 cm 104 ampoules



Small red drugs bag suitable to be fixed on a belt
Made of tear-resistant and water-repellent material
Comfortably worn on a belt with a heavy-duty belt loop
2 parts with adjustable elastics to hold 14 different sized vials securely in place
Central padding adds extra protection for fragile vials
Dimensions: about 13x12x6 cm

BOR89336	SMALL DRUGS BAG WITH ELASTICS, FOR BELT
----------	---

"Easy Rail" and "Rail" ampoules holder systems, made of painted metal with internal rubber, are available in different dimensions: visit our website www.boscarol.it for the complete list.



Book-opening drugs case empty
Made in tearproof impermeable red fabric material
2 rear Velcro-strips for the fixing in bags and emergency rucksacks
2 removable plates with elastics suitable for 28 different sized ampoules
Central padding for ampoules record
Dimensions: about 25x12x7(h) cm

BOR89330 DRUGS BAG WITH ELASTICS



Book-opening drugs case with EASY RAIL system
Made of tearproof impermeable Tyrrpol material
EASY RAIL system allows a quick insertion and removal of phials
Suitable for 60 different sized ampoules in addition to little vials (48 easy rail rigid spaces and 12 elastics)
Internal padded divider
Transparent frontal window
Dimensions: about 29x25x8 cm

BOR89340 EASY RAIL SYSTEM DRUGS CASE



Book-opening drugs case
2 removable plates able to contain 104 ampoules
8 RAIL ampoules 13 positions holders
Interior padded divider with pockets for needles and further medical material
Made in tearproof impermeable material
2 Velcro-strips for the fixing in bags or emergency rucksacks
Dimensions: about 46x8x33h cm
Weight: about 2 kg

BOR89320 MAXI DRUGS CASE RAIL



Big book-opening drugs bag
2 parts with adjustable elastics to feature 60 different sized ampoules
Central padding with transparent pocket for ampoules list and expiry date
Made in red, tear-resistant, water-resistant material
2 external Velcro-strips for the fixing in bags or rucksacks
Dimensions: about 25x25x7h cm

BOR89332 BIG DRUGS BAG WITH ELASTICS